

Product Catalog

NanoGenerator™ Nanoparticle Synthesis





NanoGenerator™ Nanoparticle Synthesis System

Nanoparticle synthesis by microfluidic technology presents advantages over the conventional batch synthesis processes due to its superior control of size and shape.

PreciGenome's NanoGenerator™ applies microfluidic approaches to synthesize nanoparticles in a continuous mode. The systems have been widely used in various applications in the drug delivery field, such as lipid nanoparticles (LNP), liposomes, PLGA nanoparticles etc.

Microfluidic Mixing System

- Controllable particle size
- Low PDI
- High encapsulation efficiency
- High reproducibility

Payloads

- DNA/mRNA/siRNA
- Proteins and peptides
- Small molecule drugs
- Other payloads

System Benefits

High Performance & Efficiency



- Tunable size (40-200nm)
- Low PDI
- High encapsulation efficiency

Open Platform



- Reagents
- Microfluidic chips

Scalable Throughput



- Flex S: 0.2-2 mL
- Flex M: 1-12 mL
- Pro/Max: 200mL/1L+

Simple Operation



- Easy setup
- Compact size
- Intuitive UI w/ touchscreen

Cost Effective



- Affordable configuration
- Low cost consumables

Custom design & OEM



- Research collaboration
- Custom design
- OEM & Contract manufacturing

| Catalog # | Name |
|-----------|---------------------------|
| PG-SYN-FS | NanoGenerator™ Flex-S |
| PG-SYN-FM | NanoGenerator™ Flex-M |
| PG-SYN-F | NanoGenerator™ Flex (S+M) |
| PG-SYN-P | NanoGenerator™ Pro |
| PG-SYN-G | NanoGenerator™ Max (GMP) |

Email: USSales@precigenome.com
Tel: +1-408-708-4602

Address: 2176 Ringwood Ave., San Jose, CA, USA
Visit us at www.precigenome.com

NanoGenerator™ System Models

Discovery & Screen



NanoGenerator™ Flex

Scientific research, screening and discovery for quick low-volume preparation

Flex-S: 0.2 to 2 mL throughput

Flex-M: 1 to 12 mL throughput



Screen & Develop




NanoGenerator™ Pro

Preclinical Studies and Development

Pro: 2 to 200 ml throughput

Scale-up



NanoGenerator™ Max (cGMP version)

Clinical development
cGMP certified manufacturing

Max: 1L throughput

Max+: >10 L throughput
(available soon)

Custom solution



NanoGenerator™ OEM

Custom design and OEM solutions
GMP certified manufacturing

High volume: >20 L throughput

High throughput: 32/64 samples per run

System Applications

Nucleic Acid Lipid Nanoparticles

- mRNA vaccines
- Rare genetic diseases
- Gene & cell therapy
- CAR-T therapeutics

Liposomes

- Cancer therapy
- Vaccine Adjuvant
- Antimicrobial therapy
- Cosmetics

Polymer Nanoparticles

- Cancer chemotherapy
- Immunology & vaccines
- Insulin delivery for diabetes

NanoGenerator™ Flex-S

Nanoparticle Synthesis System




NanoGenerator™ Flex-S nanoparticle synthesis system

PreciGenome's NanoGenerator™ Flex-S design is for small scale production. The throughput volume range is from 0.2 to 2 mL, which is perfect for formulation screening and early discovery applications.

Smaller output volume (<0.2ml per run) is achievable by recipe optimization.

| Model | NanoGenerator™ Flex-S |
|--------------------------|--|
| Mixing Cartridge | CHP-MIX-4 |
| Throughput | 0.2 to 2 mL |
| Total Flow Rate | 3 mL/min |
| Flow Rate Ratio (W:O) | 3:1 |
| Size Range | 40 to 200 nm |
| PDI | 0.05 to 0.2 |
| Encapsulation Efficiency | 85-95% |
| Payloads | DNA, mRNA, siRNA, protein, small molecules |

eGFP mRNA LNP Delivery to Jurkat Cells



Jurkat Cells transfected with Formulation #9. Green fluorescence image at 48 hours post transfection.

Example of Formulation Screening by Flex-S

| Screening Panel | | | LNP Characterization | | | Cell Study |
|-----------------|-----------------|-----------|----------------------|-------|-----|----------------|
| Formulation | Ionizable Lipid | N/P Ratio | Size (nm) | PDI | EE% | GFP expression |
| #1 | 40% | 3.57 | 56.6 | 0.19 | 86% | + |
| #2 | 40% | 5.35 | 79.9 | 0.246 | 84% | + |
| #3 | 40% | 8 | 75.2 | 0.214 | 85% | ++ |
| #4 | 60% | 5.35 | 128.5 | 0.13 | 81% | NA |
| #5 | 40% | 5.35 | 62.8 | 0.186 | 90% | ++ |
| #6 | 40% | 8 | 54.3 | 0.184 | 93% | ++ |
| #7 | 50% | 8 | 79 | 0.155 | 88% | + |
| #8 | 50% | 11 | 82.2 | 0.126 | 90% | NA |
| #9 | 50% | 8 | 87.5 | 0.12 | 91% | +++ |

Email: USSales@precigenome.com

Tel: +1-408-708-4602

Address: 2176 Ringwood Ave., San Jose, CA, USA

Visit us at www.precigenome.com/nanoparticle-synthesis

NanoGenerator™ Flex-M

Nanoparticle Synthesis System




NanoGenerator™ Flex-M nanoparticle synthesis system


PreciGenome's NanoGenerator™ applies microfluidic approaches to synthesize nanoparticles (LNP, liposome, PLGA, etc.) in a continuous mode. NanoGenerator™ Flex-M system provides a wide throughput range from 1 to 12 mL, which meets a variety of applications from early screening to animal study.

Flex-M system also provides in-line dilution option to reduce ethanol concentration instantly. It further stabilizes LNP products.


| Model | NanoGenerator™ Flex-M |
|-----------------------------|-----------------------|
| Mixing Cartridge | CHP-MIX-4 |
| Throughput | 1 to 12 mL |
| Total Flow Rate | 1 to 3 mL/min |
| Flow Rate Ratio (W:O) | 2:1 to 5:1 |
| In-line Dilution (optional) | 0.5:1 to 2:1 |
| Size Range | 40 to 200 nm |
| PDI | 0.05 to 0.2 |
| Encapsulation Efficiency | 85-95% |



mRNA LNP Synthesis



HepG2 cells were successfully transfected by eGFP mRNA LNP. Cell nucleuses were stained with Hoechst 33342 dye (blue color) before imaging.



Email: USSales@precigenome.com
 Tel: +1-408-708-4602
 Address: 2176 Ringwood Ave., San Jose, CA, USA
 Visit us at www.precigenome.com/nanoparticle-synthesis

Case Study:

T Cell Transfection by mRNA Lipid Nanoparticles

Since the first FDA approval of chimeric antigen receptor (CAR) T cell therapy in 2017, T cell engineering is continuously the hottest research field in immunotherapy and cell therapy. Current CAR T cell engineering methods use viral transductions, which induce permanent CAR expression and have potential safety concerns. To overcome this concerns, researchers are highly interested in non-viral gene delivery methods.



Recently, CAR mRNA lipid nanoparticles (LNPs) in T cell engineering have been widely studied. The transient transduction feature of mRNA LNP make it a safer profile than viral vectors. The size, homogeneousness and mRNA encapsulation efficiency are the key factors for efficient T cell transfection. Using PreciGenome's NanoGenerator™ system, customer can produce mRNA LNPs with well controlled size, high homogeneousness and excellent encapsulation efficiency.

The following data shows the size and PDI of GFP mRNA lipid nanoparticles synthesized by NanoGenerator™ Flex system. The transfection efficiency to K562 and HepG2 cell lines and human primary T cells are presented in Figure 2 and 3.




Figure 1. GFP-LNP Synthesized by PreciGenome's NanoGenerator™ Flex-S. Average sizes is 67.3 nm. PDI is 0.106.





Figure 2. GFP expression in K562 (left) and HepG2 (right) cell lines after 48 hours treated by GFP-LNP synthesized by PreciGenome's NanoGenerator™ Flex-S.

Figure 3. GFP(+) positive population of control (non-treat) and eGFP LNP treated primary T cells at 16, 40 and 64 hours. Cells were stained (1:50) using Biogend 7-AAD Viability Staining for 10 minutes. Gating: First select for individual cells (excluding doublets). Then select for the healthy cell population. Then select for viable cells by excluding cells which are positive for 7-AAD. Gate for FitC-A channel (GFP)






NanoGenerator™ Pro


Nanoparticle Synthesis System

PreciGenome's NanoGenerator™ Pro is an integrated nanoparticle synthesis system. NanoGenerator™ Pro contains more powerful pressure control modules, which provides higher throughput from 2 to 200 mL. The total flow rate range is from 4 to 12 mL/min. The flow rate ratio (W:O) is between 2:1 to 5:1.


Tunable Size




Batch Consistency




Nucleic Acids LNP Synthesis



High Encapsulation Efficiency



Cell Transfection using GFP mRNA LNP



Fluorescence Field

Bright Field


| Model | NanoGenerator™ Pro |
|--------------------------|--|
| Mixing Cartridge | CHP-MIX-3 |
| Throughput | 2 to 200 mL |
| Total Flow Rate | 4-12 mL/min |
| Flow Rate Ratio (W:O) | 2:1 to 5:1 |
| Size Range | 40 to 200 nm |
| PDI | 0.05 to 0.2 |
| Encapsulation Efficiency | Up to 99% |
| Payloads | DNA, mRNA, siRNA, protein, small molecules |

NanoGenerator™ Max


Nanoparticle Synthesis System (cGMP Version)

PreciGenome's NanoGenerator™ applies microfluidic approaches to synthesize nanoparticles in a continuous mode. The system has been widely used in the drug delivery field for a variety of applications, such as synthesis of lipid nanoparticles (LNP), liposomes, PLGA, etc.

NanoGenerator™ Max (cGMP version) is designed for clinical and commercial production (100mL-1L, 10L+). With PreciGenome's microfluidics technology, customer can easily and seamlessly transfer their early discovery results (NanoGenerator™ Flex, Pro) to late stage production (NanoGenerator™ GMP).




Cell Transfection using GFP DNA LNP



Fluorescence Field

Bright Field

High Encapsulation Efficiency



| Model | NanoGenerator™ Max | NanoGenerator™ Max+ |
|--------------------------|--|--|
| Throughput | 100mL to 1L | >= 10L |
| Total Flow Rate | >1L/hour | >10L/hour |
| Size Range | 40 to 200 nm | 40 to 200 nm |
| PDI | 0.05 to 0.2 | 0.05 to 0.2 |
| Encapsulation Efficiency | Up to 99% | Up to 99% |
| Payloads | DNA, mRNA, siRNA, protein, small molecules | DNA, mRNA, siRNA, protein, small molecules |
| Inline Dilution | Optional | Optional |

Email: USSales@precigenome.com
Tel: +1-408-708-4602

Address: 2176 Ringwood Ave., San Jose, CA, USA
Visit us at www.precigenome.com/nanoparticle-synthesis

LipidFlex™

Flexible Lipid Nanoparticle Formulation


LipidFlex™ is a 3-component lipid nanoparticle formulation that compatible with various cationic/ionizable lipids for nucleic acid encapsulation and cell transfection. LipidFlex™ Pack kit includes ionizable lipid (SM102).

- Flexible cationic/ionizable lipid ratio
- Flexible with various N/P ratio
- High nucleic acid encapsulation efficiency
- High mammalian cell transfection rate




| Model | LipidFlex™ | LipidFlex™ Pack |
|-----------------|---|--|
| Catalog # | PG-SYN-LF1ML | PG-SYN-LF1MLP |
| Components | Structural Lipid/ Cholesterol/Stabilizer | SM102/Structural Lipid/Cholesterol/ Stabilizer |
| Product size | 1000 µL | 1000 µL |
| LipidFlex Conc. | 30 mM | 30 mM |
| Ionizable lipid | NA | SM102 (20mg) |

Size & PDI




Encapsulation Efficiency




LipidFlex™ Pack Kit Experiment: HepG2 Cell Transfection Efficiency

Sample




DNA LNP, PreciGenome NanoGenerator
SM102/PG-LipidFlex = 40/60 mol%

Positive control



Lipofectamine™ 3000 (Thermo Fisher)

Negative control




Non-treat

LipidFlex™ T Cell Kit

High Efficient mRNA LNP Formulation for T Cell Transfection

LipidFlex™ T cell kit is a high efficient lipid formulation to synthesize mRNA lipid nanoparticles for primary human T cell gene delivery. Using NanoGenerator™ Flex-S system and CHIP-MIX-4 cartridge, customers can prepare potent mRNA LNP in a convenient and efficient way.



- Narrow size distribution of mRNA LNP
- High transfection efficiency
- High protein expression level
- High cell viability
- Time efficient synthesis process

| Component | Size | Storage |
|----------------------------|--------|----------|
| LipidFlex T Lipid mix | 200 µL | -20 °C |
| Formulation Buffer 1 (10x) | 60 µL | 4 - 8 °C |
| Formulation Buffer 2 | 1 mL | 4 - 8 °C |

Customer Service

- **Formulation design (Lipid NP, Liposome or PLGA)**

Customize liposomes design based on our clients' demand by varying lipid compositions, vesicle sizes, surface charges, etc.

- **Payload encapsulation**

Customize protocols to encapsulate drug molecules into lipid nanoparticle or PLGA with high encapsulation efficiency.

- **Cell study**

Cell in vitro transfection service.

- **Analysis and characterization**

Run comprehensive analysis assays for liposomes before and after encapsulation, which includes visual appearance, size distribution, stability, entrapment efficiency, encapsulation efficiency analysis, in vitro release profile analysis, release rate, etc.

Ordering Information

| Platforms | | |
|--------------------------------------|--|-----------------|
| NanoGenerator™ Flex-S | NanoGenerator™ system, throughput: 0.2mL-1mL | PG-SYN-FS |
| NanoGenerator™ Flex-M | NanoGenerator™ system, throughput: 1mL-12mL | PG-SYN-FM |
| NanoGenerator™ Flex | NanoGenerator™ system, throughput: 0.2mL-12mL, S and M module | PG-SYN-F |
| NanoGenerator™ Pro | NanoGenerator™ system, throughput: 2mL-50mL, integrated instrument | PG-SYN-P |
| NanoGenerator™ Max | GMP version nanoparticle synthesis system (available soon) | PG-SYN-G |
| Microfluidic Cartridge & Consumables | | |
| Mixing cartridge | Microfluidic mixing chip for Flex, 4 devices per chip | CHP-MIX-4 |
| Mixing cartridge | Microfluidic mixing chip for Pro, 3 devices per chip | CHP-MIX-3 |
| Mixing cartridge | Microfluidic glass mixing chip for Flex-M, 1 device per chip | CHP-MIX-G1 |
| Reservoir Connectors | Reagent Reservoirs for Flex-S, 20sets/pack | PG-MRC-SYNS-Q20 |
| Flex-S Gasket | Gaskets for Flex-S, 20pcs/pk | PG-GSK-SYNS-Q20 |
| O-ring Gasket | O-ring gasket for NanoGenerator™ Pro, 50pcs/pk | PG-ORN-SYNP-Q50 |
| Reagents (optional) | | |
| LipidFlex™ kit | 3 components lipids mixture, 1mL | PG-SYN-LF1ML |
| LipidFlex™ Pack kit | Pack of LipidFlex™ and ionizable lipid (SM102) | PG-SYN-LF1MLP |
| LipidFlex™ T cell kit | T cell transfection kit | PG-SYN-LFT |
| LipidDemo | Included in the package of the instrument | PG-SYN-LFD |
| Accessories & Service (optional) | | |
| Flex-S flow unit | Flow unit for Flex-S | PG-SYN-MNTS |
| Flex-M flow unit | Flow unit for Flex-M | PG-SYN-MNTM |
| Tubing and connector kit | Standard tuning and connectors for Flex-M | KIT-TUB-FIT-FM |
| Inline dilution kit | Inline dilution device, tubing and connectors for Flex-M | KIT-INL-DIL-FM |
| Extended Warranty | 1 to 3 years. | PG-WTY-1Y |



PreciGenome is located in the heart of Silicon Valley, San Jose, California, USA. We have been focusing on developing nanoparticle synthesis systems and solutions for our customers since we started our business. Our technology enables rapid prototyping with high quality and reliable performance for lipid nanoparticles, liposomes, PLGA, etc.

HEADQUARTER

PreciGenome LLC

Tel: 1-408-708-4602

Email: USSales@precigenome.com

Address: 2176 Ringwood Ave. San Jose, California, USA

TAIWAN DISTRIBUTOR

Ding Fong Scientific Co., Ltd

Tel: 886-4-23595717

Email: ding-fong@ding-fong.com.tw

CHINA DISTRIBUTOR

Suzhou Rainsure Biotech Co. Ltd

Tel: 86-512-67503398

Email: info@rainsurebio.com


KOREA DISTRIBUTOR

InSung Chroma-Tech Co., Ltd.

Tel: 82-226441991

Email: KRSales@precigenome.com

Some of Our Customers



Visit us at www.precigenome.com/nanoparticle-synthesis