



# NanoGenerator® MAX Nanoparticle Synthesis System




PreciGenome

Jan 2025



# What are Lipid Nanoparticles?



Lipid nanoparticles (LNPs) are self-assembling structures of natural or synthetic lipids in an aqueous environment.



# RNA-LNP Therapeutics and Vaccines




Name	Disease	Encoded antigen	Administration route	ClinicalTrials.gov identifier	Phase
<b>Infections</b>					
mRNA-1273	SARS-CoV-2	Spike	i.m.	NCT04470427	III (EUA and CMA)
BNT162b2	SARS-CoV-2	Spike	i.m.	NCT04368728	III (EUA and CMA)
CVnCoV	SARS-CoV-2	Spike	i.m.	NCT04652102	III
mRNA-1647	Cytomegalovirus	Pentameric complex and B glycoprotein	i.m.	NCT04232280	II
mRNA-1388	Chikungunya virus	Chikungunya virus antigens	i.m.	NCT03325075	I
CV7202	Rabies virus	G glycoprotein	i.m.	NCT03713086	I
<b>Cancer</b>					
mRNA-5671/V941	Non-small-cell lung cancer, colorectal cancer, pancreatic adenocarcinoma	KRAS antigens	i.m.	NCT03948763	I
mRNA-4157	Melanoma	Personalized neoantigens	i.m.	NCT03897881	II
mRNA-4650	Gastrointestinal cancer	Personalized neoantigens	i.m.	NCT03480152	I/II
HARE-40	HPV-positive cancers	HPV oncoproteins E6 and E7	i.d.	NCT03418480	I/II

Kiaie, S.H., Majidi Zolbanin, N., Ahmadi, A. *et al.* Recent advances in mRNA-LNP therapeutics: immunological and pharmacological aspects. *J Nanobiotechnol* **20**, 276 (2022).


[Nature Reviews Materials](#) volume 6, pages1078–1094 (2021)

# Lipid Nanoparticle Synthesis Methods



## Conventional Methods


### A Film hydration



- Established method
- Versatile method

- High consumption of the organic solvent
- High PDI
- Lack of reproducibility
- Need for additional downsizing step
- Difficulties in scaling-up


### B Solvent injection



- Simple and fast
- Scaling-up possibility
- Controllable

- Exposing to organic solvent
- High PDI
- Stability problem

### C Extrusion




- Established method
- Uniform and homogenous formulation

- Possible clogging of the membrane pores
- Difficulties in scaling up

Nanomaterials, Volume 11, 2021, 3440

## Microfluidic Mixing




- Low PDI
- Fast & Scalable
- Reproducibility
- Controllable
- Low cost

- Complex fabrication
- Potential clogging

VS.

# NanoGenerator® - Nanoparticle Synthesis System



## FLEX-S



## FLEX-S PLUS



- Flex-S: 0.1 – 2 ml
- Flex-S Plus: High-throughput discovery & screening

## FLEX-M




- Flex-M: 1 – 12 ml

## FLEX-M PREMIUM



- Flex-M Premium: 2– 200 ml


## MAX



Clinical development  
GMP certified manufacturing

- MAX: 50ml – 1L
- MAX (40L/h): >20L

## OEM



Custom design and OEM solutions  
GMP certified manufacturing

- >400 L throughput

DISCOVERY & SCREEN


PRE-CLINICAL DEVELOPMENT

CLINICAL  
DEVELOPMENT

CUSTOM SOLUTION




# NanoGenerator® MAX — Intro



- The NanoGenerator® Max is designed for clinical and commercial production. Two versions are available:
  - RUO: Preclinical applications
  - cGMP: Clinical and commercial production
- Two flow kits are available with different supported throughput:
  - 4.8 L/h flow kit: 50 mL – 1 L
  - 40 L/h flow kit: >20 L




# NanoGenerator® MAX — Spec



Model	NanoGenerator® MAX			
	RUO flow kit 4.8 L/h	GMP flow kit 4.8 L/h	RUO flow kit 40 L/h	GMP flow kit 40 L/h
cGMP compliance	N/A	Yes	N/A	Yes
Software (21 CFR Part 11 compliant)	Optional	Yes	Optional	Yes
Throughput	50 ml – 1 L		> 20 L	
Total flow rate	1.2 – 4.8 L/h		Up to 40 L/h	
Flow rate ratio	1:1 – 9:1		1:1 – 5:1	
Inline dilution	1:1 – 5:1			
Size range	40 – 200 nm			
PDI	0.05 – 0.2			
Encapsulation efficiency	Up to 99%			
Payload	DNA, mRNA, siRNA, protein, small molecules, etc.			
Dimension (L × W × H)	620 × 380 × 430 mm			
Weight	50 Kg		65 Kg	

# NanoGenerator® MAX — Contents



## Instrument:


- Pneumatic system
- Valves
- Flow rate sensors
- Consumable kit
- Monitor (optional)
- Pumps (optional)

## Consumable Kit:

(Sterilized, Nuclease free, pre-assembled)




- Sample bottle (aqueous)
- Sample bottle (solvent)
- Sample bottle (dilution)
- Waste bottle
- Bioprocessing bag (collection)
- Tubing & connectors
- Mixing chip






# NanoGenerator® MAX — Software



## Software (21 CFR Part 11) Features:

- Experimental parameter setting
- Experimental recipe save/load
- Real-time pressure/flow rate chart
- Historic experimental parameter tracking
- Historic pressure/flow rate tracking
- System self-diagnostic system
- Real-time flow rate diagnostic system
- Warning system
- Manual & automatic emergency stop system
- User management
- Audit trail
- Zero flow calibration
- Flow sensor maintenance & re-calibration (Service)




# NanoGenerator® MAX — Software



The screenshot displays the 'LNP™ System' configuration window. It includes fields for 'Name' (default recipe), 'Tank Size' (Water 250ml, Ethanol 250ml), 'Total flow rate' (48.0 ml/min), 'Flow rate ratio' (Aqueous: 3, Solvent: 1), 'Total volume' (4.0 ml), 'Dilution Factor' (2), and 'remove initial reagent' (3.2 ml). On the right, there are visual representations of reservoirs for Aqueous (50 ml), Solvent (50 ml), and Aqueous (3.0 ml), along with a 'Total Volume' of 4.0 ml. The 'Estimated Total Time' is 40 s. A 'Save' button is at the bottom right.

Easy-to-use UI to set parameters including:

- Total flow rate
- Flow rate ratio
- Production volume
- Inline dilution factor
- Waste volume




Easy-to-use real-time flow rate /pressure chart including:


- Flow rate
- Pressure
- Air flow rate

All parameters are tracked for aqueous, solvent, and inline dilution lines

# NanoGenerator® MAX — Performance




Nanoparticle Size vs. Total Flow Rate



- Nanoparticle size decreases as total flow rate increases
- Size decrease experiences diminishing returns when the flow rate reaches 48 ml/min

Reagents	
Aqueous phase	Sodium acetate buffer (100mM, pH5.2)
Solvent phase	LipidFlex, 15mM in ethanol


Nanoparticle Size uniformity (50 ml/fraction)



- Throughout the entire production run, there is no significant difference in the nanoparticle size and PDI

Reagents	
Aqueous phase	Phosphate-Buffered Saline (1X, pH7.4)
Solvent phase	LipidDemo, 15mM in ethanol


# NanoGenerator® — Scale Up



DISCOVERY &  
SCREEN

PRE-CLINICAL  
DEVELOPMENT

CLINICAL  
DEVELOPMENT






- Nanoparticle size is consistent across different production volumes if using optimal flow rates
- Mixing mechanism is the same for all PreciGenome instruments
- Production can be scaled up from discovery & screening to preclinical & clinical trial production

Reagents	
Aqueous phase	Sodium acetate buffer (100mM, pH5.2)
Solvent phase	LipidFlex, 15mM in ethanol




# NanoGenerator® — Scale Up



Reagents	
Aqueous phase	Sodium acetate buffer (100mM, pH5.2)
Payload	RNA (~600 nt)
Solvent phase	LipidFlex RNA-LNP kit


# Case Study: mRNA LNPs for T cell Transfection




## eGFP mRNA Lipid Nanoparticles

Z-Average Diameter: 67.3 nm

PDI: 0.106




**Figure 1.** mRNA(eGFP)-LNP Synthesized by NanoGenerator. Average diameter is 67.3 nm. PDI is 0.106. Encapsulation efficiency is 94.5% (Ribo Green RNA Quantification Kit).




**Figure 2.** GFP(+) positive population of control (non-treat) and EGFP mRNA LNP treated primary T cells at 16, 40 and 64 hours. Cells were stained (1:50) using Biolegend 7-AAD Viability Staining for 10 minutes. Gating: First select for individual cells (excluding doublets). Then select for the healthy cell population. Then select for viable cells by excluding cells which are positive for 7-AAD. Gate for FitC-A channel (GFP)


# Case Study: Bi-specific Antibody Delivered by mRNA-LNP



# Case Study: Bi-specific Antibody Delivered by mRNA-LNP



## Treatment of Two Late-Stage Breast Cancer Cases




PET-CT 12/14/2023  
Pre-treatment

PET-CT 6/18/2024  
Post-treatment

Unresectable breast cancer with skin lesion

### Case 1 Three photos

1. Appearance before treatment
2. Considerable change in appearance on skin lesions after first treatment
3. Continued improvement on skin lesions after two treatments



4/1/2024 Baseline

4/22/2024 After one Injection

Triple negative breast cancer


### Case 2

**Left:** CT scan showing a stage 3 invasive ductal carcinoma that did not respond to prior immunotherapy


**Right:** After one treatment, the tumor has dramatically resolved.



# Case Study: Bi-specific Antibody Delivered by mRNA-LNP



## Cholangiocarcinoma with Liver Metastasis



02/28/2024 after one injection




07/24/2024, after four injections

**Case Information:** A 45-year-old male patient, HBsAg positive for over 2 years, presented with intrahepatic lesions and abdominal distension. A recent CT scan revealed a large abnormal density in the liver's right lobe, enlarged abdominal lymph nodes, and a portal vein defect, indicating hepatocellular carcinoma with lymph node metastasis and portal vein cancer thrombosis. The patient's liver function was Child-Pugh grade A with some blood count abnormalities.

**The posttreatment CT scan showed dramatic shrink of the intrahepatic cholangiocarcinoma** after four intratumoral injections of the EpCAM-CD3-Fc+IM-1+IM-2 cocktail mRNA-LNP. The Patient requested more injections on 9/12/2024

# Case Study: Bi-specific Antibody Delivered by mRNA-LNP




## Liver Metastases from Colorectal Cancer

Lesion 1



**37.0 x 28.0 mm**


3/25/2024



**27.0 x 24.0 mm**


4/27/2024

Lesion 2



**70.0 x 60.0 mm**

3/25/2024




**56.0 x 43.0 mm**

4/27/2024

The enhanced CT scan of the upper abdomen showed that the intrahepatic tumor had shrunk

# Why PreciGenome?



## High Performance & Efficiency



- Tunable size (40-200 nm)
- Low PDI (0.05-0.2)
- High encapsulation efficiency

## Open Platform



- Upgradable system
- Transferable microfluidic chips

## Scalable Throughput



- Low volume for screening (Flex-S)
- Medium volume production (Flex-M)
- High volume production (Pro, MAX-GMP)

## Simple Operation



- Simple setup
- Compact size
- Intuitive UI w/ touchscreen

## Cost Effective




- Affordable configuration
- Lower cost per run

## Custom Support



- Demo, Training and Support
- Extended Warranty
- Hot swap option
- Local US company






- **Manual**
- **Standard Operation Procedure (SOP)**
- **Warranty (1 year)**
- **Documentation related to cGMP compliance (cGMP version)**
  - ✓ Installation qualification, operational qualification, performance qualification
  - ✓ Report of consumable items
  - ✓ Chemical compatibility report of consumable items
  - ✓ Report of endotoxin test
  - ✓ Report of RNase/DNase free test
  - ✓ Report of sterilization test
  - ✓ Report of ethylene oxide residue test
  - ✓ 21 CFR Part 11 report
  - ✓ Electromagnetic compatibility report
  - ✓ Report of safety regulations
  - ✓ Other reports by requesting



# Appendix II



PurePower Medical Suzhou Purepower Medical Technology Co., Ltd.

### Sterility Test Report

No.: QT/QG/Q01-01 No.: 20230724-01

Product Name	Type	Lot No.
GMP Consumable Bag	Sample	Sample

Sample ID	Sterilization Lot No.	Test Date
1, 2, 3, 4	2023072401	2023-07-24

Inoculation Method: ☒ Membrane Filter ☐ Direct Inoculation

FTM Lot No.: FTM-230724

TSB Lot No.: TSB-230724

Positive Strain: Staphylococcus aureus

Test Result:

Culture Medium	Sample ID	1	2
FTM 30-35°C	1	-	-
	2	-	-
	3	-	-
	4	-	-
TSB 20-25°C	1	-	-
	2	-	-
	3	-	-
	4	-	-

Incubation Temperature (°C): 37

Dish No.: 1

Incubation Time: 24h, 48h

Average:

Conclusion: ☒ Con

Remark:

Notice: Put the "P" in the "C", in the "R"

Tested by/Date: Yn Yan 2023

PurePower Medical Suzhou Purepower Medical Technology Co., Ltd.

### EO Residual Test Report

No.: QT/QG/Q04-02 No.: 20230724-01

Sample Name	Lot No.	Sample No.
GMP Consumable	Sample	Sample

Test Date: 2023-07-24

Test Reference: JL-QG/Q04 EO Residual Test

Item	Sample1	Sample2	Sample3	Conclusion	Remark

Tested by/Date: Yn Yan 2023

CTI PHARMA

### Extractables Test Report

Report title: The Extractables Study Report of Nanoparticle Synthesis System Consumables Kit

Report number: EL-REP-23-019.01-E

Project No.: N/A

Customer: PreciGenome LLC

Address: 2176 Ringwood Ave. San Jose, CA, 95131, USA

Testing laboratory: Centre Testing International Pinchuang (Shanghai) Co., Ltd.

Testing laboratory address: 1351 Wanfang Road, Minhang District, Shanghai

Underwriters Laboratories (UL LLC) IEC/EN Safety Report

Model: PG-SYN-G

Device Description: NanoGenerator™ Max Nanoparticle Synthesis System

Applicant: PreciGenome LLC

Manufacturer: Same as Applicant

Manufacturing Facility(ies): Suzhou PreciGen Unit 202, Building 2, Suzhou, 215125

Report No.: E526160-D1003

Report (Re) Issue Date: 2023-12-06

Base Standard(s): EN 61010-1:2011

Additional Standards: N/A

Report Types: This report contains information for the safety evaluation of the product.

This report covers the Safety evaluation above.

Test Report issued under the responsibility of:

UL Solutions

### TEST REPORT IEC 61010-1

Safety requirements for electrical equipment for measurement, control, and protection

Report Number: E526160-D1003

Date of issue: 2023-12-06

Total number of pages: 1

Name of Testing Laboratory preparing the Report: UL Solutions

Applicant's name: PreciGenome LLC

Address: 2176 Ringwood Ave, San Jose, CA, 95131, USA

Test specification: IEC 61010-1

Standard: IEC 61010-1

Test procedure: In accordance with the standard

Non-standard test method: None

TRF template used: IEC 61010-1

Test Report Form No.: IEC 61010-1

Test Report Form Originator: UL Solutions

Master TRF: IEC 61010-1

Copyright © 2021 IEC System of Conformity Assessment Standards (IECEE System). All rights reserved.

This publication may be reproduced in whole or in part without prior approval of the IEC System of Conformity Assessment Standards (IECEE) provided that the reader's interpretation of the reproduced material is the responsibility of the manufacturer, after completion of an EC Declaration of Conformity and compliance with all relevant EC Directives.

This report is not valid as a CB Test Report and appended to a CB Test Certificate.

General disclaimer: The test results presented in this report are for information only. This report shall not be reproduced, except in full, without the prior approval of the General Manager. This verification is subjected to the governance of the General Conditions of Services, printed overleaf.

UL Solutions

### Statement of Compliance

UL-CCIC Company Limited  
No. 2, Chengwan Road, Suzhou Industrial Park, Suzhou 215122, China  
T: +86-512-6808 6400  
F: +86-512-6808 4099

Project No.: 4790895205-2.1

Applicant: PreciGenome LLC

Address of Applicant: 2176 Ringwood Ave, San Jose, CA, 95131, USA

Product Description: NanoGenerator™ Max Nanoparticle Synthesis System

Model No.: PG-SYN-G

Test Standard: EN IEC 61326-1:2021

Test Report Number(s): 4790895205-2.1-1

This verification of EMC Compliance has been granted to the applicant based on the results of the tests, performed by laboratory of UL-CCIC Company Ltd. on the sample of the above-mentioned product in accordance with the provisions of the relevant specific standards and Directive 2014/53/EU. The CE mark as shown below can be used, under the responsibility of the manufacturer, after completion of an EC Declaration of Conformity and compliance with all relevant EC Directives.

CE

Leon Wu  
Leon Wu  
UL-CCIC Company Limited.

Copyright of this verification is owned by UL-CCIC Company Ltd. and may not be reproduced other than in full and with the prior approval of the General Manager. This verification is subjected to the governance of the General Conditions of Services, printed overleaf.

